

**Vivinex™ Toric**

CLARITY OF VISION AND  
OUTSTANDING ROTATIONAL STABILITY



**Vivinex™ Toric preloaded** in the **Vivinex™ iSert® injector**  
provides **outstanding delivery consistency**

**HOYA**  
SURGICAL OPTICS

# Vivinex™ Toric

## Clarity of vision and outstanding rotational stability


Designed for outstanding optical quality, Vivinex™ Toric has proven rotational stability for precise astigmatism correction and provides patients with an astigmatic cornea with clarity of vision. Product quality, dedication and attention to detail are deeply rooted in our Japanese heritage, and with 2 million lenses implanted worldwide, surgeons' trust in Vivinex™ is proven.

### Vivinex™ Toric:

- Glistening-free hydrophobic acrylic IOL material<sup>[1,2]</sup>
- Proprietary aspheric optic design for improved image quality<sup>[3]</sup>
- Active oxygen processing treatment, a smooth surface and square optic edge to reduce PCO<sup>[1,4,5,6,7,8,9,10]</sup>
- Median rotation 1.1° (range 0.0° – 5.0°)  
100% of lenses (n=103) had ≤ 5° of rotation from their initial axis at end of surgery through all follow up visits at 1 hour, 1 week, 1 month and 6 months<sup>[11]</sup>

### Textured-rough haptic surface for better grip

Vivinex™ haptics have different surface structures. They provide better grip inside the capsular bag and are designed to reduce the potential for adhesion to the optic surface.



Rough haptic surface

Textured haptic edge


Rough haptic surface

# Reliable outcomes through outstanding rotational stability<sup>[11]</sup>


Early clinical results with the Vivinex™ preloaded IOL platform at the University of Vienna confirm outstanding rotational stability. In a group of 103 eyes, no lens rotated more than 5 degrees from orientation at the end of surgery to 6 months postoperatively.

## The Vivinex™ IOL platform shows outstanding rotational stability between surgery and one week post op, without outliers, when compared to AcrySof\*, Tecnis\* and enVista\*.<sup>[11,12]</sup>

Maximum rotation values observed in the first week following surgery<sup>[11,12]</sup>



In the timeframe of up to 6 months post-op, no implanted Vivinex™ IOL rotated more than 5° from initial axis, in comparison to 4.8% of the AcrySof\* IOL, 7.5% of the Tecnis\* IOL and 13.9% of the enVista\* IOL.<sup>[11,12]</sup>




Most studies published about rotational stability of toric IOLs do not measure early rotation from the initial axis at end of surgery. The baseline for rotation measurements is often the axis of alignment at 1 day postoperative.<sup>[11,12]</sup>

## Proprietary aspheric optic design for improved image quality


In the presence of decentration Vivinex™ minimises coma when compared with other aspheric IOLs at 4.0 mm pupil diameter.<sup>[3]</sup>

Studies have shown that the mean decentration of an IOL following cataract surgery is  $0.4 \pm 0.2$  mm with a range up to 1.7 mm.<sup>[13]</sup>

### Reduced coma caused by off-axis alignment





### Reduced impact on optical quality caused by off-axis alignment<sup>[3]</sup>



# More than 9 million HOYA preloaded IOLs implanted worldwide

Vivinex™ iSert® provides outstanding delivery consistency and offers the following features:<sup>[14,15]</sup>

- Injector tip outer diameter of 1.70 mm
- Uni-directional slider
- Screw injector with fixed grip



## Step A

Infuse the OVD into the injector through the infusion port. Fill up the area indicated by dotted lines.



## Step B

Press the release tabs, lift up and remove the cover from the case.



## Step C

Hold body with thumb and push the slider slowly forward until it stops. Remove the injector from the case.




## Step D

Carefully insert the injector tip into the eye through the incision, keeping the slit of the tip in a downward position. Slowly rotate the injector knob clockwise, to inject the lens into the capsular bag.

The handling shown above illustrates in summary the product application and does not replace the Instruction For Use.

# Introducing the HOYA Toric Calculator 4



According to several studies, standard toric IOL calculations tend to result in a calculation error of astigmatic value. Most often in:<sup>[16]</sup>

- undercorrection of against-the-rule astigmatism (ATR)
- overcorrection of with-the-rule astigmatism (WTR).

The HOYA Toric Calculator (HTC) 4 can take account of posterior corneal astigmatism in the calculation by giving the option to apply the Abulafia-Koch Regression formula.


The Abulafia-Koch Regression, applied to a clinical patient cohort, has been shown to improve predictability of TIOIOL refractive outcomes.<sup>[16]</sup>

For Vivinex™ Toric IOL calculation please visit [www.HOYAtoric.com](http://www.HOYAtoric.com)

## The HTC 4 at a glance

- Choice of three different cylinder power options allows the doctor to select the most suitable IOL model based on residual astigmatism and axis
- Optional Abulafia-Koch Regression formula can account for the added astigmatic effect of the posterior cornea when measured by standard keratometry of the anterior corneal surface
- Option to display calculation results as plus (+) or minus (-) cylinder
- Adjustable keratometer index (default 1.3375)
- Optional axial length data entry to adjust the cylinder power of the toric IOL at the corneal plane
- Numerous different print and export options are available with customized orientation of the printed eye image

# Technical characteristics



Vivinex™ Toric	
Model name	XY1A
Optic design	Biconvex with square, thin and textured optic edge Anterior: Aspheric design Posterior: Toric design
Optic & haptic materials	Hydrophobic acrylic Vivinex™ with UV- and blue light filter
Haptic design	Textured-rough haptic surface
Diameter (optic/OAL)	6.00 mm / 13.00 mm
Power	+10.00 to +30.00 D (in 0.50 D increments)
Cylinder power <sup>[17]</sup>	1.00 to 6.00 D (T2 to T9) T2 to T3 in 0.50 D increments T3 to T9 in 0.75 D increments
Nominal A-constant**	118.9
Optimized constants***	Haigis $a_0 = -0.8028$ $a_1 = 0.2133$ $a_2 = 0.2245$
	Hoffer Q     pACD = 5.697
	Holladay 1     sf = 1.934
	SRK/T     A = 119.198
Injector	Vivinex™ iSert® preloaded
Front injector tip outer diameter	1.70 mm
Recommended incision size	2.20 mm

Model XY1A	Cylinder power at IOL plane	Cylinder power at corneal plane <sup>[18]</sup>
T2	1.00 D	0.69 D
T3	1.50 D	1.04 D
T4	2.25 D	1.56 D
T5	3.00 D	2.08 D
T6	3.75 D	2.60 D
T7	4.50 D	3.12 D
T8	5.25 D	3.64 D
T9	6.00 D	4.17 D

\*\* The A-constant is presented as a starting point for the lens power calculation. When calculating the exact lens power, it is recommended that calculations be performed individually, based on the equipment used and operating surgeon's own experience.

\*\*\* These optimized constants for the calculation of intraocular lens power published by IOLCon on their website: <https://iolcon.org> are calculated from 1,475 clinical results for Vivinex™ Model XY1/XC1 as of September 24, 2021. These constants are based on actual surgical data and are provided by IOLCon as a starting point for individual constant optimizations. The information available on the website is based on data originating from other users and not by HOYA Surgical Optics ("HSO"). HSO therefore does not warrant the correctness, completeness and currentness of the contents on the said website.

- 1 HOYA data on file. DoF-CTM-21-002, HOYA Medical Singapore Pte. Ltd, 2021
  - 2 Tandogan, T. et al. (2021): In-vitro glistening formation in six different foldable hydrophobic intraocular lenses. In BMC Ophthalmol 21, 126
  - 3 Pérez-Merino, P.; Marcos, S. (2018): Effect of intraocular lens decentration on image quality tested in a custom model eye. In: Journal of cataract and refractive surgery 44 (7), p. 889–896.
  - 4 Leydolt, C. et al. (2020): Posterior capsule opacification with two hydrophobic acrylic intraocular lenses: 3-year results of a randomized trial. In: American journal of ophthalmology 217 (9), p. 224–231.
  - 5 Giacinto, C. et al. (2019): Surface properties of commercially available hydrophobic acrylic intraocular lenses: Comparative study. In: Journal of cataract and refractive surgery 45 (9), p. 1330–1334.
  - 6 Werner, L. et al. (2019): Evaluation of clarity characteristics in a new hydrophobic acrylic IOL in comparison to commercially available IOLs. In: Journal of cataract and refractive surgery 45 (10), p. 1490–1497.
  - 7 Nanavaty, M. et al. (2019): Edge profile of commercially available square-edged intraocular lenses: Part 2. In: Journal of cataract and refractive surgery 45 (6), p. 847–853.
  - 8 Matsushima, H. et al. (2006): Active oxygen processing for acrylic intraocular lenses to prevent posterior capsule opacification. In: Journal of cataract and refractive surgery 32 (6), p. 1035–1040.
  - 9 Farukhi, A. et al. (2015): Evaluation of uveal and capsule biocompatibility of a single-piece hydrophobic acrylic intraocular lens with ultraviolet-ozone treatment on the posterior surface. In: Journal of cataract and refractive surgery 41 (5), p. 1081–1087.
  - 10 Eldred, J. et al. (2019): An In Vitro Human Lens Capsular Bag Model Adopting a Graded Culture Regime to Assess Putative Impact of IOLs on PCO Formation. In: Investigative ophthalmology & visual science 60 (1), p. 113–122.
  - 11 Schartmüller, D. et al. (2019): True rotational stability of a single-piece hydrophobic intraocular lens. In: The British journal of ophthalmology 103 (2), p. 186–190.
  - 12 Schartmüller, D. et al. (2020): Comparison of Long-Term Rotational Stability of Three Commonly Implanted Intraocular Lenses. In American journal of ophthalmology 220, pp. 72–81.
  - 13 Harrer, A. et al. (2017): Variability in angle  $\kappa$  and its influence on higher-order aberrations in pseudophakic eyes. In: Journal of cataract and refractive surgery 43 (8), p. 1015–1019.
  - 14 Oshika, T.; Wolfe, P. (2019): In vitro comparison of delivery performance of 4 preloaded intraocular lens injector systems for corneal and sclerocorneal incisions. In: Journal of cataract and refractive surgery 45 (6), p. 840–846.
  - 15 Haldipurkar, S. et al. (2020): Incision size changes after cataract surgery with intraocular lens implantation: comparison of 2 preloaded IOL implantation injectors. In: Journal of cataract and refractive surgery 46 (2), p. 222–227.
  - 16 Abulafia, A. et al. (2016): New regression formula for toric intraocular lens calculations. In: Journal of cataract and refractive surgery 42 (5), p. 663–671.
  - 17 At IOL plane.
  - 18 Based on an average pseudophakic human eye.
- \* Third-party trademarks used herein are the property of their respective owners.

Information contained is intended for health care professionals. For a full list of indications and contra indications please refer to the Instructions For Use. Some of the products and/or specific features as well as the procedures featured in this document may not be approved in your country and thus may not be available there. Design and specifications are subject to change without prior notice as a result of ongoing technical development. Please contact our regional representative regarding individual availability in your country. HOYA, Vivinex and iSert are trademarks of the HOYA Corporation or its affiliates.

©2021 HOYA Medical Singapore Pte. Ltd. All rights reserved.



HOYA Medical Singapore Pte. Ltd | 455A Jalan Ahmad Ibrahim | Singapore 639939



HOYA Surgical Optics GmbH | De-Saint-Exupéry-Straße 10 | 60549 Frankfurt am Main | Germany  
Hotline DE: Tel. +49 (0)800 664 2 664 | Fax +49 (0)800 774 2 774

[hoyasurgicaloptics.com](http://hoyasurgicaloptics.com)



2021-09-27\_HSOE\_XY1A\_BR\_EN

**HOYA**  
SURGICAL OPTICS